

RB-1000

- 10.4" Full Touch Screen
- Active & Passive Trainer with the adjustable recumbent seat
- Available to swivel the seat by 45°.
- Display and Setting
 - Speed
 - Resistance
 - Time
 - Direction
- Spasm control



Model Name	SP-RB-1000
Code Number	100081
Description	Active & Passive Trainer for Lower Body with the adjustable recumbent seat
Dimension of Device (L x W x H, mm)	1600 x 530 x 1210
Weight	78kg
Power Supply	110/220V 50-60Hz
Power Consumption	96VA
Monitor	10.4" LCD Touch Screen
Dimension of Monitor	264 x 220 x 39
Training	Leg



SP-1100 | SP-2100 | SP-3100

SP Standard Series

- Active & Passive Trainer
- Voice Guidance
- Display and Setting
 - Speed
 - Resistance
 - Time
 - Direction
- Height Adjustable
- Spasm control

| SP-1100



| SP-2100



| SP-3100



Model Name	SP-1100	SP-2100	SP-3100
Code Number	100004	100005	100006
Description	Motor Fit Upper & Lower Body with 7.1" Touch Screen	Motor Fit Lower Body with 7.1" Touch Screen	Motor Fit Upper Body with 7.1" Touch Screen
Dimension of Device (L x W x H, mm)	920 x 550 x 1080	730 x 520 x 1070	920 x 550 x 1080
Weight	40kg	33kg	32kg
Power Supply	110/220V 50-60Hz	110/220V 50-60Hz	110/220V 50-60Hz
Power Consumption	96VA	96VA	96VA
Monitor	7" LCD Touch Screen	7" LCD Touch Screen	7" LCD Touch Screen
Dimension of Monitor	200 x 227 x 55	200 x 227 x 55	200 x 227 x 55
Training	Arm & Leg	Leg	Arm